




ALTAMIRA INSTRUMENTS

The First Name in Custom Reactor Systems

BenchCAT™ X000R HP Hybrid

A CUSTOM CATALYST TESTING INSTRUMENT



- Performs both reaction testing *AND* catalyst characterization
- Multiple workstations (2, 3, 4, 5, 6, or 8 workstations)
- Fully automated workstations
- Independent experiments on each workstation
- Designed for unattended operation

Reaction Testing:

- Temperature Programmed Reactions
- Isothermal reaction studies
- Pressures from ambient up to 50 bar
- Furnace temperatures up to 1200°C
- Analysis via auxiliary MS or GC
- Quartz, 316 SS, or custom reactor tube
- Sample holder sizes to 25 mL

Catalyst Characterization:

- TPR, TPO, TPD, TPRx, Pulse Chemisorption
- Dynamic BET
- Furnace temperatures up to 1200°C
- Ambient pressure operation
- Subambient testing to -100°C
- Analysis via integral TCD or auxiliary MS (or GC)
- Quartz sample holder (to 1 gram)
- Monolith holder available (1"Ø x 2" long)

The Altamira Advantage:

- Altamira will customize this instrument to meet your exact research needs today.
- Altamira will customize this instrument in the future to meet your changing research needs.

149 Delta Drive, Pittsburgh, PA 15238 • 412-963-6385 (Voice) • 412-963-6485 (Fax)
Email: info@altamirainstruments.com

Visit our Website: www.altamirainstruments.com

Technical Specification–BenchCAT™ 4000R HP Hybrid

System Capabilities:

Catalyst Characterization and Reaction Testing: The BenchCAT 4000R HP Hybrid is a four station automated reactor system designed to:

1. Conduct temperature-programmed or isothermal experiments with catalysts loaded in a straight tube reactor at ambient or pressure conditions. External analysis will be performed with a mass spectrometer.
2. Perform characterization experiments (TPR, TPO, TPD, Pulse Chemisorption) using catalysts at ambient pressures. Analysis will be made internally, using the TCD, or externally with a mass spectrometer.
3. The instrument is configured to perform either reaction experiments or characterization experiments on an automated basis. Reaction and characterization will always be performed as separate experiments.

Basic System Specifications

A. Reaction Testing

| | |
|---------------------------------------|---|
| Number of Stations | 4 |
| System Operating Pressure | Ambient – 20 barg (higher pressures available) |
| MFC Inlet Pressure | 350 psig (Varies according to operating pressure) |
| Maximum Furnace Operating Temperature | 700°C |
| Maximum Post Reactor Line Temperature | 150°C |
| Number of MFCs (Reaction Side) | 3 (per station) |
| Number of MFCs (Shared) | 3 (per station, blended) |
| Number of Vaporized Liquid Feeds | 1 (HPLC pump per station) |
| Materials of Construction | |
| • Upstream Plumbing | 316SS |
| • Downstream Plumbing | 316SS |
| • Sample tubes | Glass lined stainless steel |
| • Sample tube option | 316 SS monolith holder |
| • Wetted Parts | 316SS |
| • Seal Materials | Premium |
| • Catalyst Charge | To 1 gram |
| Overall System Footprint | 70" x 60" x 30" (H x L x W, approximate) |
| Utility Requirements: | |
| • Oil-free, Dry Air | Multiple 80 psig utility connections |
| • System power | 208 VAC, 90A |
| • Computer power | 110V/15A |

B. Catalyst Characterization

| | |
|---------------------------------------|-----------------|
| Number of Stations | 4 |
| System Operating Pressure | Ambient |
| Gas Inlet Pressure Range | 50 psig |
| Gas Outlet Pressure Range | near ambient |
| Maximum Furnace Operating Temperature | 1200°C |
| Maximum TCD Operating Temperature | 200°C |
| Materials of Construction | |
| • Plumbing | 316SS |
| • Sample U-tubes | Quartz |
| • Wetted Parts | 316SS |
| • Seal Materials | Premium seals |
| • Catalyst Charge | 0.1 – 1 g |
| Number of Treatment/Carrier Ports | 4 (per station) |
| Number of Blend Ports | 2 (per station) |

